




**Double stand holder**  
For NEUROWERK EEG/EMG combined systems.



**Keyboard & mouse disinfected**  
For use on the ICU.



**Height-adjustable mounting for TFT monitor**



**Slots for swivel arms/stands**  
For attachment of EEG/EMG headbox, photic stimulator, etc.



With laptop

**Different castors and castor sizes**  
With different qualities, e.g. conductive



**With additional metal baskets/shelves**




**With ergonomic sliding handle**



**With additional shelves for printers**



**With additional drawers/leaves**



**With separate mouse leaf**



## EQUIPMENT OPTIONS

## MEDICAL TROLLEY AND EQUIPMENT




manufactured by  
**SIGMA Medizin-Technik GmbH**

Business Park Am Gruendel 2  
D-09423 Gelenau – Germany

Tel.: +49 (0) 37297 – 825 – 0  
Fax: +49 (0) 37297 – 825 – 23

www.neuowerk.de  
info@neuowerk.de



manufactured by  
**SIGMA Medizin-Technik GmbH**

Business Park Am Gruendel 2  
D-09423 Gelenau – Germany

Tel.: +49 (0) 37297 – 825 – 0  
Fax: +49 (0) 37297 – 825 – 23

www.neuowerk.de  
info@neuowerk.de



## Type H The Flexible One

- » Metal/plastic trolley with 4 movable castors and locking brake (Castors: 100mm and 75mm)
- » Trolley dimensions: 550mm (W), 650mm (D), 1000mm (H)\*
- » With sliding handle, storage area and built-in socket strip, integrated connection for potential compensation cable
- » MBT holder for swivel arm or stand
- » With bracket for isolating transformer

### OPTIONS:

- » With NEUROWERK EEG and/or EMG equipment
- » With desktop PC or mini PC with TFT monitor (rotating and swivelling) or with laptop
- » Expandable with additional shelves, metal baskets (front/back) as well as various accessories
- » Individually and flexibly configurable



## Type M The Ergonomic One

- » Metal/plastic trolley with 4 movable castors and locking brake (Castors: 100mm or 125mm)
- » »Trolley dimensions: 610mm (W), 620mm (D), 884-1332mm (H)\*, height-adjustable (useable standing and sitting)
- » With sliding handle, storage are and built-in socket strip, integrated connection for potential compensation cable
- » MBT holder for swivel arm or stand
- » With built-in isolation transformer

### OPTIONS:

- » With NEUROWERK EEG or EMG equipment
- » With desktop PC or mini PC with TFT monitor (rotating and swivelling)
- » Expandable with additional shelves, metal baskets (front/back) as well as various accessories
- » Individually and flexibly configurable

## Type E The Solid One

- » Space-saving wooden/plastic trolley with 4 movable castors and locking brake (Castors: 75mm or 100mm)
- » Trolley dimensions: 508mm (W), 540mm (D), starting from 780mm (H) \*
- » With integrated isolating transformer and trolley keyboard leaf
- » With back for cable cover

### OPTIONS:

- » With NEUROWERK EEG and/or EMG equipment
- » With swivel arm (EMG) or mounted stand (EEG)
- » With desktop PC and TFT monitor or with laptop
- » Expandable with printer or drawer and side metal basket and various accessories
- » Individually and flexibly configurable



## Type I The Robust One

- » Metal/plastic trolley with 4 movable castors and locking brake (Castors: 100mm or 125mm)
- » Trolley dimensions: 555mm (W), 625mm (D), from 840mm (H)\*
- » With integrated isolating transformer

### OPTIONS:

- » With NEUROWERK EEG or EMG equipment
- » With swivel arm (EMG) or mounted stand (EEG)
- » With desktop PC and TFT monitor (rotating and swivelling) or with laptop
- » Expandable with additional shelves and drawers and various accessories
- » Individually and flexibly configurable

